

Preparation of Pancreas Lysate

1. Pancreas organs were removed from the animal and flash frozen in liquid nitrogen. Organs were kept at 80°C until use.
2. The frozen tissue was suspended in 5 volumes of ice-cold lysis buffer containing 50mM Tris (pH 7.4), 5mM EDTA (pH 8), 1% Triton X-100, and protease inhibitor cocktail (Complete, EDTA-free, Roche).
3. The tissue was immediately homogenized with a Polytron homogenizer. Homogenates were then centrifuged at 150,000 *g* for 30 minutes at 4°C .
4. The resulting supernatant was transferred to a clean tube and protein content was established. Protein concentration was adjusted to 3 mg/ml with lysis buffer. The lysates were kept at -80°C until use.